

ixi Lite B7

The first step into smart care



Product Overview

Elevate the standard of care with ixi Lite. From indoor alarms to accurate location tracking, this device is designed to provide safety, peace of mind, and a deeper connection between caregivers and residents. Requires ixi Gateways.

Key Features

- **Indoor Alarm Button:** Ensure quick response in emergencies with the indoor alarm button. Instantly notify caregivers or responders when assistance is needed, enhancing safety within enclosed spaces.
- **Indoor Location Tracking with Beacons:** With our innovative beacon technology, ixi Lite offers precise indoor location tracking. Always be informed about the whereabouts of residents, even in complex indoor environments.
- **Real-Time Resident Monitoring:** Gain insights into resident activity patterns and behaviour. Monitor for any abnormalities or changes, enabling proactive caregiving and timely intervention.
- **Waterproof Design:** Designed to withstand everyday activities, ixi Lite is waterproof, providing uninterrupted safety even in wet conditions.
- **Extended Battery Life:** The device boasts a battery life of over half a year, minimising the need for frequent recharging and ensuring continuous safety coverage.
- **Versatile Wear Options:** Wearable as a wristband or pendant, ixi Lite adapts to individual preferences, making it comfortable for residents to carry at all times.
- **Automatic System Updates:** Stay up to date effortlessly. ixi Lite's system updates occur

automatically, eliminating the need for manual intervention and ensuring optimal functionality.

Technical specification

General specifications

Dimensions (H x W x D) 35,9 x 35,9 x 11mm

Dimensions wristband (H x W x D) 267 x 15 x 1 mm

Weight approx: 22g

Waterproof: IP67

Battery: CR2032 (220mAh)

Radios

Frequency hopping 2.4Ghertz radio technology

Monitoring features

Advertise iBeacon & Eddystone & 3-axis sensor data

Detachable wristband

Monitoring and control

ixicare installer dashboard

Firmware upgrade

Device configuration

Package content

ixi Lite B7

Wristband

Device manual

User interfaces

Push-button

Regulatory compliance

iBeacon MFi License (iBC-14-00582)

FCC Regulation (FCC Part 15.)

CE Regulations (Included EN300328/301489/60950/62479)